

The Wave

A Bathing & Transfer System

Hygiene and health go together. When we make bathing easier, our clients are healthier. The Wave system has been thoroughly reimagined to improve bathing and showering for the client and for the caregiver.

Responding to requests from therapists and caregivers, we widened the seat to accommodate larger clients. And to protect skin integrity, we added generous padding under soft knit fabric.

In addition to the tub stand and rolling shower stand, we've added an ingenious new tub transfer base to completely change the way you transfer and bathe large clients. There's nothing like it.



Find letters of medical necessity and informative articles at: www.rifton.com/wave

The Wave bath chair is terrific!
In the past, lifting a child into
a tub has been a back killer,
but the new transfer base
eliminates this hazard entirely.
The casters allowing lateral
movement in the tub are a
welcome addition, and the fact
that I can easily collapse the
base and hang it on the wall is
a fabulous space saver.

Linda K. Miller, DPT Oueensbury, New York



Unique Features & Benefits

Include these in your letters of medical necessity:

- **Extra wide seat** with open sides, 15" width for small and medium, 17½" for the large.
- Large size accommodates adults up to 6'2" and 250 lbs.
- **Completely configurable**; all sizes of the Rifton Wave fit on any of the 3 base options.
- Optional calf rest adjusts independently of seat angle.
- **Tub transfer base** enables users—after being placed in the bath chair using a mechanical lift—to be rolled back and then rotated over the tub for showering.
- Conversion kit allows user to economically move up from small to medium size by just purchasing a new back, seat and cover.
- **Least costly alternative** to remodeling a bath room or assisted living, without resorting to bed baths.

Features of the Wave

No tools required



Generously sized

The wide seat and open design enables large users to fit comfortably. Seat width on the large size is 17½".



One-hand adjustment

Adjust backrest angle with one hand, keeping the other hand on the client for safety and reassurance.



Full leg support

The optional calf rest provides full support for the lower legs. Adjust it independently of the seat by pushing the gray buttons and rotating up or down. It must be installed at the time of purchase.






Typical bath tub height range

——18 in.

———14 in.







Required components

Bath chair size and color



Choice of chest strap



Without lateral positioning

The simple chest strap safely positions the user in the Wave bath chair. It attaches with hook and loop where needed along the full height of the back, and buckles securely behind.





With lateral positioning

Rifton designed this chest strap with adjustable lateral straps to give additional lateral support helping users maintain better trunk alignment in the bath chair. The soft fabric laterals will not get in the way during transfers.





Accessories



Head blocks

Optional head blocks provide lateral support and positioning for the user's head.





Leg straps

Use the fully adjustable leg straps to position your client's hips in abduction (or just to provide additional security).





Calf rest

The new calf rest option provides more support for the lower legs and is especially needed when the bath chair is used on one of the stands.

Note: Calf rest must be chosen at the time of the original purchase. It cannot be added later.





Tub stand

The optional tub stand raises the bath chair either $5\frac{1}{2}$ " or $11\frac{1}{2}$ " higher in the tub depending on the position of the tub stand legs.





Shower stand

The shower stand lets you roll your client into the shower stall at a height that allows for better cleaning. Better cleaning means better health. Locking casters (4") provide smooth, easy maneuverability.



Accessories (continued)

Tub transfer base

Bathing equipment for people with disabilities tends to be big and bulky,

taking up far too much space in bathrooms shared by the whole family. Now using Rifton's ingenious new tub transfer base you can perform healthy transfers without crowding your bathroom with permanently mounted equipment. When bath time is over the base collapses and hangs compactly on the wall.



Unique two-stage movement



The tub transfer base starts in the position farthest out of the tub. It is then rolled back over the tub, and the brakes are set.



As the upper portion with the bath chair is rotated 90° it moves approximately 12° back over the tub in one motion.



The tub transfer base fits bath tubs up to 24" high by adjusting the leg length to clear the edge of the tub.





1. With the client secured in the chair, push the transfer base to the back of the tub as far as it will go and lock the wheels.



2. Press either of the gray pivot latches and rotate the bath chair and client over the tub.



3. As the bath chair rotates, it continues to move back over the tub. Gently lift the client's legs to clear the edge of the tub.



4. In the final position the bath chair is over the tub and can be adjusted as desired for showering.

Free up the tub.

Rifton's folding, storable tub transfer base makes it possible for everyone in the family to still use the tub.

When you are finished, remove the base from the tub, fold the legs and hang it on the wall hook provided.





Dimensions of the Wave

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40-56	50-74



Key user dimension: height

Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth.

Important: Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.

Item dimensions (inches)	Z210 small	Z220 medium	Z230 large				
Seat depth	10½	13	17				
Seat width	15	15	171/2				
Back height	23	26½	32				
Seat height	2–7	2–7	2–7				
Seat height (incl. tub stand)	12½-18½	121/2 – 181/2	121/2 - 181/2				
Seat height (incl. shower stand)	271/2	271/2	271/2				
Seat height (incl. tub transfer base) with long leg inserts	21-26 24½-29½	21-26 24½-29½	21-26 24½-29½				
Length of mesh area when flat (without calf rest)	35	41	50½				
Length of mesh area when flat (with calf rest)	421/2	49	61				
Width of base	18	18	21				
Length of base when flat	291/2	291/2	291/2				
Tub stand height (legs flat)	8	8	8				
Tub stand length (legs flat)	39	39	39				
Tub stand width (legs flat)	17½	171/2	171/2				
Shower stand height	22	22	22				
Shower stand width	27	27	27				
Shower stand length	34	34	34				
Tub transfer base footprint	28 x 30	28 x 30	28 x 30				
Max tub height for transfer base with long leg inserts	20 24	20 24	20 24				
Maximum working load (lbs)	75	150	250				
Back and seat angle adjustments	ack and seat angle adjustments 0°, 15°, 30°, 45°, 60°, 75°						



Without calf rest (Z269)

The conversion kit

(two styles shown above) is the smart solution for growing children. Convert the small to a medium bath chair for a fraction of the cost of a new chair.

With calf rest (Z267)

We'll send you a new larger fabric cover, back and seat frames and chest strap with instructions on how to remove the old and install the new parts on your existing base.



The Wave has a washable fabric cover that is easily removed for cleaning.







All three sizes are available in blue and pink.



How to buy

- **1. Measure** your client to determine the size you need. (See dimension chart for help.)
- **2. Select** the color, appropriate bath chair and chest strap.
- **3. Choose** the accessories you need.
- **4. Add** up the total cost.
- **5. Order** online in minutes at **Rifton.com** or see the back of this catalog for other ways to order.
- **6. Call** us any time at **800.571.8198** if you need assistance. We're always here to help.

The tool-free ${\bf conversion}$ ${\bf kit}$ converts your small chair to medium as your child grows.

Fabric without calf rest

Blue Z268 \$310 Pink Z269 \$310 (Includes seat, back, cover and chest strap)

Fabric with calf rest

Blue	Z266	\$385
Pink	Z267	\$385
(Includes seat, back,	cover and che	est strap)



			Z210 sr	nall	Z220 n	nedium	Z230 I	arge
Rec	Chair		Z211	\$260	Z221	\$270	Z231	\$330 🗌
Required components	Fabric cover	Blue	Z244	\$185 🗌	Z246	\$195	Z248	\$205 🗌
		Pink	Z245	\$185 🗌	Z247	\$195	Z249	\$205 🗌
pone	Chest strap	With lateral positioning	Z212	\$62 🗌	Z222	\$72 🗌	Z232	\$82 🗌
ents		Without lateral positioning	Z213	\$41 🗌	Z223	\$52 🗌	Z233	\$62 □
Opt	Head blocks		Z214	\$88 🗌	Z214	\$88	Z234	\$88
ional	Leg straps			\$41 🗌	Z225	\$52 🗌	Z235	\$62 🗌
Optional accessories	Calf rest (calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216	\$115 🗌	Z226	\$130 🗌	Z236	\$145
ssori	Tub stand		Z217	\$260 🗌	Z217	\$260	Z217	\$260
S	Shower stand		Z218	\$700 🗌	Z218	\$700	Z218	\$700 🗌
	Tub transfer base (calf rest recommended)	Counterclockwise into tub	Z227	\$825 🗌	Z227	\$825 🗌	Z227	\$825 🗌
		Clockwise into tub	Z228	\$825 🗌	Z228	\$825 🗌	Z228	\$825 🗌
	Legs for tub transfer base	Standard	Z224	\$52 🗌	Z224	\$52 🗌	Z224	\$52 🗌
	CEGS FOR LIAB GRANE BASE	Long (tub edge over 20" high)	Z238	\$52 🗌	Z238	\$52 🗌	Z238	\$52 🗌
	Additional head blocks		Z214	\$88 🗌	Z214	\$88 🗌	Z234	\$88 🗌
	Additional fabric cover		Z219	\$185 🗌	Z229	\$195	Z239	\$205 🗌
	Additional chest strap	With lateral positioning	Z212	\$62 □	Z222	\$72 🗌	Z232	\$82 🗌
		Without lateral positioning	Z213	\$41 🔲	Z223	\$52 🗌	Z233	\$62 □
Tota	I	1	Z210	\$	Z220	\$	Z230	\$